



How to Effectively Increase Capacity of Absorbent Protective Products

WHAT WORKS ...



- Combining a booster pad with a disposable brief or pull-on is a smart way to provide longer wear time and extra leakage protection while traveling, taking a long walk in the park, or for everyday use. Its flow-through design allows absorption into the booster pad, then excess fluid passes through into the garment.

WHAT DOESN'T ...



- Moisture-proof barriers in personal care pads, liners, guards and feminine hygiene pads are designed to prevent fluid from flowing through to the bottom garment. This style of protection is ideal for regular underwear, but does not increase the effectiveness of a host incontinence product.



What is a Booster Pad?

Booster Pads, also known as “booster inserts” or simply “boosters,” are products placed inside of another incontinence product

FLOW-THROUGH DESIGN

(typically disposable briefs or pull-on underwear) to provide additional absorbent capacity. Boosters feature a flow-through design that allows them to absorb an initial amount of fluid, before additional fluid passes through the product, into the host incontinence garment below it. Because of this, **they should never be used in regular underwear**, or accidents will likely occur.

Combining a booster pad with a disposable brief or pull-on is a smart way to provide longer wear time and extra leakage protection while traveling, taking a long walk in the park, or for everyday use. Booster pads can be folded and even doubled for added protection, but should never be cut.


TIP: Positioning of booster pads is very important!

Women: Place booster in the middle of host garment

Men: Place booster toward the front of host garment

No part of the booster should stick out, over the edge of the host garment’s pad, or leakage could occur.

BOOSTER PADS




VS.



What is a Personal Care Pad?

Many people confuse **booster** pads with **personal care** pads (also known as guards, liners, feminine hygiene pads or “period pads”). Tranquility offers a line of personal care pads designed for incontinence protection. These types of pads differ from boosters in that **they have a moisture-proof backing** that is designed to stop fluid from leaking through. For this reason, personal care pads are meant to be worn in regular underwear and **should not be worn inside of another incontinence product**. Doing this will likely cause excess fluid to spill off the side of the personal care pad and leak out of the host garment.

PERSONAL CARE PADS



Remember:

- **Booster Pads** – Used **only** with **incontinence garments**
- **Personal Care Pads** – Used **only** with **regular underwear**



Have questions or want a free booster sample?
1-866-865-6101 TranquilityProducts.com