

Tranquility® Essential Booster Pads

Flow-through pad designed to add capacity to any disposable incontinence product.



Previously Branded As



Features and Benefits

Flow-Through Design • Clothlike
Adhesive Strip • Latex-Free



ADDITIONAL
CAPACITY

Superabsorbent booster adds capacity to your primary product, extending your wear time.



SOFT
COMFORT

Soft and flexible design allows you to position the booster comfortably where you need it.



SECURE
PLACEMENT

Adhesive strip secures the booster to your primary product to prevent sliding, keeping you comfortable and protected.

HEAVY



(Regular)

| Item No. | Size | Length/Width | Capacity | Case Count | HCPC Code | Best Suited for Host Garment Sizes |
|----------|---------|--------------|---------------|--------------------|-----------|------------------------------------|
| 2770 | Youth | 12" x 3.25" | 5.5 oz/163 mL | 180 (6 bags of 30) | T4535 | Youth, XS |
| 2760 | Regular | 12" x 4.25" | 6.7 oz/198 mL | 200 (8 bags of 25) | T4535 | S, M, L |
| 2762 | Long | 15" x 4.25" | 8.0 oz/236 mL | 200 (8 bags of 25) | T4535 | XL, 2XL |

MODERATE



(Regular)

| Item No. | Size | Length/Width | Capacity | Case Count | HCPC Code | Best Suited for Host Garment Sizes |
|----------|---------|--------------|---------------|--------------------|-----------|------------------------------------|
| 19244 | Regular | 12" x 4" | 5.3 oz/157 mL | 200 (8 bags of 25) | T4535 | S, M, L, XL |




Want more information? 1-800-467-3224
TranquilityProducts.com/brand-evolution



What is a Booster Pad?

Booster pads have a **flow-through design** allowing them to absorb initial fluid before passing additional moisture into a primary incontinence product. Booster pads should always be placed inside another disposable incontinence product that has a moisture-proof backing. **Do not use with regular underwear.**




Did You Know?

Booster pad placement is important. **The pads can be placed in targeted areas where additional absorption is needed.**



Booster pads can be folded and even doubled for added protection, but should never be cut.



Host Garment – The pad's adhesive side should face out, away from the body. Make sure the pad is fully contained within the host incontinence product; no part should extend beyond the edges.

For Women – Place the pad more toward the bottom.

For Men – Place the pad more toward the front.



Want more information?
1-800-467-3224 • TranquilityProducts.com

