

# Usage Guidelines

## Tranquility® Personal Care Pads



### Purpose

Tranquility products are manufactured with special superabsorbent polymers. Because of their unique design, moisture is pulled into the core of the product keeping the surface against the skin dry and promoting longer wear times.

### What to expect

- Skin will stay dry because the moisture is wicked away.
- The hourglass design conforms to the body for added comfort.
- Decrease in changing clothing and linens due to leakage.
- Adhesive strips secure pad in place in regular underwear.
- Preservation of dignity for those that need some protection but do not require a disposable brief/pull-on.

### Sizes and Capacity

Item #	Description	Fluid Capacity
2380	Super 10.5" x 5.5"	400 mL
2381	Ultimate 13.5" x 6.5"	450 mL
2382	OverNight 16.5" x 7.25"	500 mL

### When to change/replace

- There is no set time between changes.
- Perform checks per facility protocol.



For additional information or questions  
Call: 1-800-467-3224 x 340 or email: [clinical@pbenet.com](mailto:clinical@pbenet.com)